


Keep track of your data over time with the Floodlight™ MS App



Floodlight MS is a science-based smartphone app that remotely keeps track of function over time to help patients and healthcare providers expand the conversation around personalized MS care. The collected data will be visible to patients through the dashboard within the mobile app. The Study Investigators will be able to see the data for their patients via the web-portal. Having access to your performance data will help you and your doctor to better understand your disease progression.



Your contribution helps shaping the future for people living with MS

Clinical studies are the reason medical care has come so far, and the smartphone monitoring app is especially important because it aims to measure the daily reality of MS. By recording your weekly progress, you're helping health professionals better understand MS and its management.

We, and the entire MS community, appreciate you taking an active role to advance MS knowledge in a way only you can. Because nobody understands your journey with MS better than you.

Floodlight MS includes **5** activities to collect data on your **hand function**, **walking ability**, and **cognition**, which are commonly affected areas of function in people living with MS.

PINCH A TOMATO **DRAW A SHAPE**

hand function

This block contains two activity icons: a tomato and a hand drawing a spiral. Below the icons are the labels "PINCH A TOMATO" and "DRAW A SHAPE". At the bottom of the block is the category label "hand function".

U-TURN TEST **2-MINUTE WALK**

walking ability

This block contains two activity icons: a U-shaped arrow and a stick figure walking. Below the icons are the labels "U-TURN TEST" and "2-MINUTE WALK". At the bottom of the block is the category label "walking ability".

MATCH THE SYMBOLS

cognition

This block contains one activity icon: two overlapping squares, one with the number 6 and one with a diagonal line. Below the icon is the label "MATCH THE SYMBOLS". At the bottom of the block is the category label "cognition".

Please see the user manuals on [FloodlightMS.com](https://www.floodlightms.com) for directions on how to properly perform the activities in the app.

The Floodlight MS app also features a Journal to record how you are feeling and any symptoms you might be experiencing. It can help guide you and your healthcare provider toward improving conversations about your care.

FL MS TONiC Study information

What is the purpose of this observational study?

The purpose of this observational study is to examine how well the FL MS App can measure any effects of MS on your symptoms and abilities, thereby improving monitoring of your MS. It will also study whether there is any relation between commonly used tests to assess disease activity in MS and data collected using the FL MS App.

The FL MS App is a science-based smartphone app designed to inform conversations and healthcare provider decisions by enabling people with MS to regularly record data on their function between clinic visits.

After completion of the run-in period, you will enter the observational period for a period of 3 to 4 years.

To participate in this study, your smartphone has to be compatible with the FL MS mobile application. Your doctor will help you to check the device compatibility during the enrollment visit.

Important when using Floodlight MS



1. **Complete the assessments** according to the defined schedule.



2. **Recharge your phone regularly** and keep it in good working condition.



3. **Inform the study center immediately in the event of damage, malfunctioning, loss or replacement of your smartphone.**



4. **Bring your personal smartphone** to every visit at the clinic.

GETTING SUPPORT

If you face any difficulties during the study, please always contact your study investigator. To protect your privacy as study participant, we do not offer direct customer support services for patients participating in clinical trials.



Roche takes patient privacy seriously. If you'd like to learn more, please refer to the Floodlight MS privacy notice.

